

产品名称: DHX8 Rabbit Monoclonal Antibody
产品货号: AMRe87397



产品概述 (Summary)

产品名称 (Production Name)	DHX8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DHX8
别名 (Alternative Names)	DDX8; Dhr2; HRH1; PRP22; PRPF22
基因 ID (Gene ID)	1659
蛋白 ID (SwissProt ID)	Q14562.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:10000
蛋白分子量 (Molecular Weight)	Calculated MW:139 kDa; Observed MW:139 kDa

研究背景 (Background)


产品名称: DHX8 Rabbit Monoclonal Antibody
产品货号: AMRe87397



This gene is a member of the DEAH box polypeptide family. The encoded protein contains the DEAH (Asp-Glu-Ala-His) motif which is characteristic of all DEAH box proteins, and is thought to function as an ATP-dependent RNA helicase that regulates the release of spliced mRNAs from spliceosomes prior to their export from the nucleus. This protein may be required for the replication of human immunodeficiency virus type 1 (HIV-1). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of DHX8 in HeLa,293T cell lysates using DHX8 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .