

产品名称: PRPF31 Rabbit Monoclonal Antibody
产品货号: AMRe87398



产品概述 (Summary)

产品名称 (Production Name)	PRPF31 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRPF31
别名 (Alternative Names)	RP11; PRP31; SNRNP61; NY-BR-99
基因 ID (Gene ID)	26121
蛋白 ID (SwissProt ID)	Q8WWY3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:100-1:200, FC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:56 kDa; Observed MW:56 kDa

研究背景 (Background)


产品名称: PRPF31 Rabbit Monoclonal Antibody
产品货号: AMRe87398




This gene encodes a component of the spliceosome complex and is one of several retinitis pigmentosa-causing genes. When the gene product is added to the spliceosome complex, activation occurs.[provided by RefSeq, Jan 2009]

研究领域 (Research Area)


图片 (Image Data)



Western blot detection of PRPF31 in C2C12, PC3 cell lysates using PRPF31 antibody (1:1000 diluted).



Immunofluorescent analysis of HeLa cells using AMRe87398 antibody (green), and DAPI (blue).



Western blot analysis of extracts from HeLa, 293T, RAW264.7, PC-12 cells using AMRe87398 at 1:3000.

产品名称: PRPF31 Rabbit Monoclonal Antibody
产品货号: AMRe87398



注意事项 (Note)

For research use only .