

产品名称: BAT5 Rabbit Monoclonal Antibody
产品货号: AMRe87413



产品概述 (Summary)

产品名称 (Production Name)	BAT5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BAT5
别名 (Alternative Names)	BAT5; NG26; PP199; D6S82E
基因 ID (Gene ID)	7920
蛋白 ID (SwissProt ID)	O95870.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:63 kDa; Observed MW:55 kDa

研究背景 (Background)


产品名称: BAT5 Rabbit Monoclonal Antibody
产品货号: AMRe87413



A cluster of genes, BAT1-BAT5, has been localized in the vicinity of the genes for tumor necrosis factor alpha and tumor necrosis factor beta. These genes are all within the human major histocompatibility complex class III region. The protein encoded by this gene is thought to be involved in some aspects of immunity. Alternatively spliced transcript variants have been described. [provided by RefSeq, Apr 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of BAT5 in HeLa, Jurkat cell lysates using BAT5 antibody (1:1000 diluted).

注意事项 (Note)

For research use only.