

产品名称: SPARCL1 Rabbit Monoclonal Antibody
产品货号: AMRe87416



产品概述 (Summary)

产品名称 (Production Name)	SPARCL1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SPARCL1
别名 (Alternative Names)	SC1; MAST9; PIG33; MAST 9
基因 ID (Gene ID)	8404
蛋白 ID (SwissProt ID)	Q14515.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:75 kDa; Observed MW:130 kDa

研究背景 (Background)


产品名称: SPARCL1 Rabbit Monoclonal Antibody
产品货号: AMRe87416



SPARCL1 (Secreted Protein, Acidic and Rich in Cysteines-like 1), also known as hevin, SC1 or MAST9, is a member of the SPARC family of extracellular glycoproteins. SPARCL1 is an anti-adhesive protein that is widely expressed in tissues such as brain, heart, lung, muscle and kidney, but not liver.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human heart tissue using AMRe87416 at 1:1000.

注意事项 (Note)

For research use only .