

产品名称: YKL-40/CHI3L1 Rabbit Monoclonal Antibody
产品货号: AMRe87417



产品概述 (Summary)

产品名称 (Production Name)	YKL-40/CHI3L1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YKL-40/CHI3L1
别名 (Alternative Names)	GP39; ASRT7; GP-39; YK-40; YKL40; CGP-39; YKL-40; YYL-40; HC-gp39; HCGP-3P; hCGP-39
基因 ID (Gene ID)	1116
蛋白 ID (SwissProt ID)	P36222.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:200,FC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:43 kDa; Observed MW:43 kDa

产品名称: YKL-40/CHI3L1 Rabbit Monoclonal Antibody
产品货号: AMRe87417



研究背景 (Background)

Chitinases catalyze the hydrolysis of chitin, which is an abundant glycopolymer found in insect exoskeletons and fungal cell walls. The glycoside hydrolase 18 family of chitinases includes eight human family members. This gene encodes a glycoprotein member of the glycosyl hydrolase 18 family. The protein lacks chitinase activity and is secreted by activated macrophages, chondrocytes, neutrophils and synovial cells. The protein is thought to play a role in the process of inflammation and tissue remodeling. [provided by RefSeq, Sep 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from THP-1 cells using AMRe87417 at 1:1000.

注意事项 (Note)

For research use only .