

产品名称: Haptoglobin Rabbit Monoclonal Antibody
产品货号: AMRe87422



产品概述 (Summary)

产品名称 (Production Name)	Haptoglobin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Haptoglobin
别名 (Alternative Names)	HP-1; preHP2
基因 ID (Gene ID)	15439
蛋白 ID (SwissProt ID)	Q61646.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:39 kDa; Observed MW:45 kDa

研究背景 (Background)


产品名称: Haptoglobin Rabbit Monoclonal Antibody
产品货号: AMRe87422



This gene encodes a plasma glycoprotein called haptoglobin that binds free hemoglobin. The encoded preproprotein undergoes proteolytic processing to generate alpha and beta subunits that form a disulfide-linked tetrameric protein that plays an important role in the sequestration and clearance of extracorporeal hemoglobin. Mice lacking the encoded protein exhibit stunted development of lymphoid organs associated with lower counts of mature T and B cells in the blood and secondary lymphoid compartments. [provided by RefSeq, Aug 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse lung tissue using AMRe87422 at 1:1000.

注意事项 (Note)

For research use only .