

产品名称: ADAM10 Rabbit Monoclonal Antibody
产品货号: AMRe87449



产品概述 (Summary)

产品名称 (Production Name)	ADAM10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADAM10
别名 (Alternative Names)	RAK; kuz; AD10; AD18; MADM; CD156c; CDw156; HsT18717
基因 ID (Gene ID)	102
蛋白 ID (SwissProt ID)	O14672.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:84 kDa; Observed MW:90,68 kDa

研究背景 (Background)


产品名称: ADAM10 Rabbit Monoclonal Antibody
产品货号: AMRe87449




Members of the ADAM family are cell surface proteins with a unique structure possessing both potential adhesion and protease domains. This gene encodes an ADAM family member that cleaves many proteins including TNF- α and E-cadherin. Alternate splicing results in multiple transcript variants encoding different proteins that may undergo similar processing. [provided by RefSeq, Feb 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from MCF-7 cells using AMRe87449 at 1:1000.



Western blot analysis of extracts from HCT116, A431 cells using AMRe87449 at 1:3000.

注意事项 (Note)

For research use only .