

产品名称: CD20 Rabbit Monoclonal Antibody
产品货号: AMRe87453



产品概述 (Summary)

产品名称 (Production Name)	CD20 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD20
别名 (Alternative Names)	B1; S7; Bp35; CD20; CVID5; MS4A2; LEU-16
基因 ID (Gene ID)	931
蛋白 ID (SwissProt ID)	P11836.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:200,ICC/IF 1:10-1:50,FC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:33 kDa; Observed MW:33 kDa

研究背景 (Background)


产品名称: CD20 Rabbit Monoclonal Antibody
产品货号: AMRe87453



This gene encodes a member of the membrane-spanning 4A gene family. Members of this nascent protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and nonlymphoid tissues. This gene encodes a B-lymphocyte surface molecule which plays a role in the development and differentiation of B-cells into plasma cells. This family member is localized to 11q12, among a cluster of family members. Alternative splicing of this gene results in two transcript variants which encode the same protein. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Raji cells using AMRe87453 at 1:1000.

注意事项 (Note)

For research use only .