

产品名称: Nuclear Matrix Protein p84 Rabbit Monoclonal Antibody
产品货号: AMRe87459



产品概述 (Summary)

产品名称 (Production Name)	Nuclear Matrix Protein p84 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Nuclear Matrix Protein p84
别名 (Alternative Names)	P84; HPR1; P84N5
基因 ID (Gene ID)	9984
蛋白 ID (SwissProt ID)	Q96FV9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:50-1:100,ICC/IF 1:100-1:200,FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:76 kDa; Observed MW:84 kDa

产品名称: Nuclear Matrix Protein p84 Rabbit Monoclonal Antibody
产品货号: AMRe87459




研究背景 (Background)


HPR1 is part of the TREX (transcription/export) complex, which includes TEX1 (MIM 606929), THO2 (MIM 300395), ALY (MIM 604171), and UAP56 (MIM 142560).[supplied by OMIM, Nov 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, 293T, HepG2 cells using AMRe87459 at 1:1000.



Western blot analysis of extracts from SW620 cells and Rat liver, Mouse brain tissue using AMRe87459 at 1:2000.

注意事项 (Note)

For research use only .