

产品名称: Integrin linked ILK Rabbit Monoclonal Antibody
产品货号: AMRe87462



产品概述 (Summary)

产品名称 (Production Name)	Integrin linked ILK Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Integrin linked ILK
别名 (Alternative Names)	P59; ILK-1; ILK-2; p59ILK; HEL-S-28
基因 ID (Gene ID)	3611
蛋白 ID (SwissProt ID)	Q13418.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000, ICC/IF 1:200-1:500, FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:51 kDa; Observed MW:51 kDa


研究背景 (Background)

产品名称: Integrin linked ILK Rabbit Monoclonal Antibody
产品货号: AMRe87462

This gene encodes a protein with a kinase-like domain and four ankyrin-like repeats. The encoded protein associates at the cell membrane with the cytoplasmic domain of beta integrins, where it regulates integrin-mediated signal transduction. Activity of this protein is important in the epithelial to mesenchymal transition, and over-expression of this gene is implicated in tumor growth and metastasis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87462 at 1:1000.

注意事项 (Note)

For research use only .