

产品名称: INMT Rabbit Monoclonal Antibody
产品货号: AMRe87464



产品概述 (Summary)

产品名称 (Production Name)	INMT Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	INMT
别名 (Alternative Names)	TEMT
基因 ID (Gene ID)	11185
蛋白 ID (SwissProt ID)	O95050.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:29 kDa; Observed MW:29 kDa

研究背景 (Background)


产品名称: INMT Rabbit Monoclonal Antibody
产品货号: AMRe87464




N-methylation of endogenous and xenobiotic compounds is a major method by which they are degraded. This gene encodes an enzyme that N-methylates indoles such as tryptamine. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the downstream FAM188B (family with sequence similarity 188, member B) gene. [provided by RefSeq, Nov 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of INMT in 3T3 cell lysates using INMT antibody(1:1000 diluted).



Western blot analysis of extracts from C6 cells using AMRe87464 at 1:1000.

注意事项 (Note)

For research use only .