

产品名称: DNAJC1 Rabbit Monoclonal Antibody
产品货号: AMRe87476



产品概述 (Summary)

产品名称 (Production Name)	DNAJC1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNAJC1
别名 (Alternative Names)	HTJ1; MTJ1; ERdj1; DNAJL1
基因 ID (Gene ID)	64215
蛋白 ID (SwissProt ID)	Q9NVH1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:64 kDa; Observed MW:64 kDa

研究背景 (Background)


产品名称: DNAJC1 Rabbit Monoclonal Antibody
产品货号: AMRe87476



The membrane protein encoded by this gene is a DNAJ-like heat shock protein that binds the molecular chaperone BiP. In addition, the encoded protein contains two SANT domains that have been shown to bind serpin alpha1-antichymotrypsin and inter-alpha trypsin inhibitor heavy chain 4. [provided by RefSeq, Jul 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87476 at 1:1000.

注意事项 (Note)

For research use only .