

产品名称: CD9 Rabbit Monoclonal Antibody
产品货号: AMRe87480



产品概述 (Summary)

产品名称 (Production Name)	CD9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD9
别名 (Alternative Names)	MIC3; MRP-1; BTCC-1; DRAP-27; TSPAN29; TSPAN-29
基因 ID (Gene ID)	928
蛋白 ID (SwissProt ID)	P21926.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:25 kDa; Observed MW:22 kDa

研究背景 (Background)

产品名称: CD9 Rabbit Monoclonal Antibody
产品货号: AMRe87480




This gene encodes a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Tetraspanins are cell surface glycoproteins with four transmembrane domains that form multimeric complexes with other cell surface proteins. The encoded protein functions in many cellular processes including differentiation, adhesion, and signal transduction, and expression of this gene plays a critical role in the suppression of cancer cell motility and metastasis.
[provided by RefSeq, Jan 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87480 at 1:1000.



Western blot analysis of extracts from U87-MG cells using AMRe87480 at 1:1000.

注意事项 (Note)

For research use only .