

产品名称: Phospho-NF-kB p65 (Ser536) Rabbit Monoclonal Antibody
产品货号: AMRe87488



产品概述 (Summary)

产品名称 (Production Name)	Phospho-NF-kB p65 (Ser536) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-NF-kB p65
别名 (Alternative Names)	p65; NFKB3
基因 ID (Gene ID)	5970, 19697, 309165
蛋白 ID (SwissProt ID)	Q04206.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:60 kDa; Observed MW:65 kDa

产品名称: Phospho-NF-kB p65 (Ser536) Rabbit Monoclonal Antibody
产品货号: AMRe87488




研究背景 (Background)

NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFKB1 or NFKB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFKB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87488 at 1:1000.

注意事项 (Note)

For research use only .