

**产品名称: BAFF Rabbit Monoclonal Antibody**  
**产品货号: AMRe87494**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | BAFF Rabbit Monoclonal Antibody        |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB, ICC/IF, FC, IP                     |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                                      |
|------------------------|----------------------------------------------------------------------|
| 基因名 (Gene Name)        | BAFF                                                                 |
| 别名 (Alternative Names) | DTL; BAFF; BLYS; CD257; TALL1; THANK; ZTNF4; TALL-1; TNLG7A; TNFSF20 |
| 基因 ID (Gene ID)        | 10673                                                                |
| 蛋白 ID (SwissProt ID)   | Q9Y275.                                                              |

## 产品应用 (Application)

|                          |                                                                   |
|--------------------------|-------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000, ICC/IF 1:100-1:200, FC 1:200-1:500, IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | Calculated MW:31 kDa; Observed MW:31 kDa                          |

## 研究背景 (Background)


**产品名称: BAFF Rabbit Monoclonal Antibody**  
**产品货号: AMRe87494**



The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptors TNFRSF13B/TACI, TNFRSF17/BCMA, and TNFRSF13C/BAFFR. This cytokine is expressed in B cell lineage cells, and acts as a potent B cell activator. It has been also shown to play an important role in the proliferation and differentiation of B cells. Alternatively spliced transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Mar 2011]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from U937 cells using AMRe87494 at 1:1000.

## 注意事项 (Note)

For research use only .