

产品名称: TAPP-1 Rabbit Monoclonal Antibody
产品货号: AMRe87497



产品概述 (Summary)

产品名称 (Production Name)	TAPP-1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAPP-1
别名 (Alternative Names)	TAPP1
基因 ID (Gene ID)	59338
蛋白 ID (SwissProt ID)	Q9HB21.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:46 kDa; Observed MW:46 kDa

研究背景 (Background)

产品名称: TAPP-1 Rabbit Monoclonal Antibody
产品货号: AMRe87497




This gene encodes a pleckstrin homology domain-containing adapter protein. The encoded protein is localized to the plasma membrane where it specifically binds phosphatidylinositol 3,4-bisphosphate. This protein may be involved in the formation of signaling complexes in the plasma membrane. Polymorphisms in this gene are associated with age-related macular degeneration. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 5.[provided by RefSeq, Sep 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of TAPP-1 in K562 cell lysates using TAPP-1 antibody(1:1000 diluted).



Western blot analysis of extracts from HeLa cells using AMRe87497 at 1:1000.

注意事项 (Note)

For research use only .