

产品名称: NEK6 Rabbit Monoclonal Antibody
产品货号: AMRe87512



产品概述 (Summary)

产品名称 (Production Name)	NEK6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEK6
别名 (Alternative Names)	SID6-1512
基因 ID (Gene ID)	10783
蛋白 ID (SwissProt ID)	Q9HC98.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:36 kDa; Observed MW:36 kDa

研究背景 (Background)


产品名称: NEK6 Rabbit Monoclonal Antibody
产品货号: AMRe87512




The protein encoded by this gene is a kinase required for progression through the metaphase portion of mitosis. Inhibition of the encoded protein can lead to apoptosis. This protein also can enhance tumorigenesis by suppressing tumor cell senescence. Several transcript variants encoding a few different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of NEK6 in HeLa, C6 cell lysates using NEK6 antibody (1:1000 diluted).



Western blot analysis of extracts from U2OS cells and Rat spleen tissue using AMRe87512 at 1:1000.

注意事项 (Note)

For research use only.