

产品名称: PP2A alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe87528



产品概述 (Summary)

产品名称 (Production Name)	PP2A alpha/beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PP2A alpha/beta
别名 (Alternative Names)	RP-C; PP2Ac; PP2CA; PP2Calpha
基因 ID (Gene ID)	5515
蛋白 ID (SwissProt ID)	P67775.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:200,ICC/IF 1:20-1:50,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:36 kDa; Observed MW:36 kDa

研究背景 (Background)


产品名称: PP2A alpha/beta Rabbit Monoclonal Antibody
产品货号: AMRe87528



This gene encodes the phosphatase 2A catalytic subunit. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. This gene encodes an alpha isoform of the catalytic subunit. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of PP2A alpha/beta in K562,Rat Brain,C6,3T3 cell lysates using PP2A alpha/beta antibody(1:1000 diluted).

Western blot analysis of extracts from HepG2, RAW264.7 cells and Rat muscle tissue using AMRe87528 at 1:1000.

注意事项 (Note)

For research use only .