

产品名称: DHX36 Rabbit Monoclonal Antibody
产品货号: AMRe87530



产品概述 (Summary)

产品名称 (Production Name)	DHX36 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DHX36
别名 (Alternative Names)	RHAU; Ddx36; mKIAA1488; 2810407E23Rik
基因 ID (Gene ID)	72162
蛋白 ID (SwissProt ID)	Q8VHK9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:114 kDa; Observed MW:114 kDa

研究背景 (Background)


产品名称: DHX36 Rabbit Monoclonal Antibody
产品货号: AMRe87530



Enables DNA helicase activity; histone deacetylase binding activity; and nucleic acid binding activity. Involved in several processes, including G-quadruplex DNA unwinding; positive regulation of cell differentiation; and regulation of transcription by RNA polymerase II. Acts upstream of or within response to exogenous dsRNA and response to virus. Located in cytosol. Is expressed in heart and telencephalon. Orthologous to human DHX36 (DEAH-box helicase 36). [provided by Alliance of Genome Resources, Apr 2022]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse brain tissue using AMRe87530 at 1:1000.

注意事项 (Note)

For research use only .