

产品名称: ETS2 Rabbit Monoclonal Antibody
产品货号: AMRe87532



产品概述 (Summary)

产品名称 (Production Name)	ETS2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ETS2
别名 (Alternative Names)	ETS2IT1
基因 ID (Gene ID)	2114
蛋白 ID (SwissProt ID)	P15036.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:53 kDa; Observed MW:53 kDa

研究背景 (Background)


产品名称: ETS2 Rabbit Monoclonal Antibody
产品货号: AMRe87532



This gene encodes a transcription factor which regulates genes involved in development and apoptosis. The encoded protein is also a protooncogene and shown to be involved in regulation of telomerase. A pseudogene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from RAW264.7 cells using AMRe87532 at 1:1000.

注意事项 (Note)

For research use only .