

产品名称: SLAMF6 Rabbit Monoclonal Antibody
产品货号: AMRe87564



产品概述 (Summary)

产品名称 (Production Name)	SLAMF6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLAMF6
别名 (Alternative Names)	KALI; NTBA; CD352; KAL1b; Ly108; NTB-A; SF2000
基因 ID (Gene ID)	114836
蛋白 ID (SwissProt ID)	Q96DU3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:500-1:1000, FC 1:200-1:500, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:37 kDa; Observed MW:60 kDa

研究背景 (Background)


产品名称: SLAMF6 Rabbit Monoclonal Antibody
产品货号: AMRe87564



The protein encoded by this gene is a type I transmembrane protein, belonging to the CD2 subfamily of the immunoglobulin superfamily. This encoded protein is expressed on Natural killer (NK), T, and B lymphocytes. It undergoes tyrosine phosphorylation and associates with the Src homology 2 domain-containing protein (SH2D1A) as well as with SH2 domain-containing phosphatases (SHPs). It functions as a coreceptor in the process of NK cell activation. It can also mediate inhibitory signals in NK cells from X-linked lymphoproliferative patients. Alternative splicing results in multiple transcript variants encoding distinct isoforms.[provided by RefSeq, May 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of SLAMF6 in Ramos cell lysates using SLAMF6 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .