

产品名称: Glycogen Synthase 1 Rabbit Monoclonal Antibody
产品货号: AMRe87580



产品概述 (Summary)

产品名称 (Production Name)	Glycogen Synthase 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Glycogen Synthase 1
别名 (Alternative Names)	GSY; GYS
基因 ID (Gene ID)	2997
蛋白 ID (SwissProt ID)	P13807.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW:84 kDa; Observed MW:90 kDa


研究背景 (Background)

产品名称: Glycogen Synthase 1 Rabbit Monoclonal Antibody
产品货号: AMRe87580


The protein encoded by this gene catalyzes the addition of glucose monomers to the growing glycogen molecule through the formation of alpha-1,4-glycoside linkages. Mutations in this gene are associated with muscle glycogen storage disease. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Sep 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Glycogen synthase 1 in HeLa cell lysates using Glycogen synthase 1 antibody(1:500 diluted).



Western blot analysis of extracts from Mouse heart tissue using AMRe87580 at 1:5000.

注意事项 (Note)

For research use only .