

产品名称: DAP kinase 3 Rabbit Monoclonal Antibody
产品货号: AMRe87582



产品概述 (Summary)

产品名称 (Production Name)	DAP kinase 3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DAP kinase 3
别名 (Alternative Names)	DLK; ZIP; ZIPK
基因 ID (Gene ID)	1613
蛋白 ID (SwissProt ID)	O43293.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:53 kDa; Observed MW:53 kDa

研究背景 (Background)


产品名称: DAP kinase 3 Rabbit Monoclonal Antibody
产品货号: AMRe87582




Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of ZIP Kinase in HeLa, A549, MCF-7, U251 cell lysates using ZIP Kinase antibody (1:1000 diluted).



Western blot analysis of extracts from T24 cells using AMRe87582 at 1:1000.

注意事项 (Note)

For research use only .