

产品名称: Staufen Rabbit Monoclonal Antibody
产品货号: AMRe87584



产品概述 (Summary)

产品名称 (Production Name)	Staufen Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Staufen
别名 (Alternative Names)	STAU; PPP1R150
基因 ID (Gene ID)	6780
蛋白 ID (SwissProt ID)	O95793.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:100-1:200, FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:63 kDa; Observed MW:63 kDa

研究背景 (Background)


产品名称: Staufen Rabbit Monoclonal Antibody
产品货号: AMRe87584




Staufen is a member of the family of double-stranded RNA (dsRNA)-binding proteins involved in the transport and/or localization of mRNAs to different subcellular compartments and/or organelles. These proteins are characterized by the presence of multiple dsRNA-binding domains which are required to bind RNAs having double-stranded secondary structures. The human homologue of staufen encoded by STAU, in addition contains a microtubule-binding domain similar to that of microtubule-associated protein 1B, and binds tubulin. The STAU gene product has been shown to be present in the cytoplasm in association with the rough endoplasmic reticulum (RER), implicating this protein in the transport of mRNA via the microtubule network to the RER, the site of translation. Five transcript variants resulting from alternative splicing of STAU gene and encoding three isoforms have been described. Three of these variants encode the same isoform, however, differ in their 5'UTR. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Staufen in K562, C6, 3T3 cell lysates using Staufen antibody (1:1000 diluted).



Western blot analysis of extracts from HepG2 cells using AMRe87584 at 1:2000.

注意事项 (Note)

For research use only.

产品名称: Staufen Rabbit Monoclonal Antibody
产品货号: AMRe87584

