

产品名称: TIAM2 Rabbit Monoclonal Antibody
产品货号: AMRe87589



产品概述 (Summary)

产品名称 (Production Name)	TIAM2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TIAM2
别名 (Alternative Names)	STEF; TIAM-2
基因 ID (Gene ID)	26230
蛋白 ID (SwissProt ID)	Q8IVF5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:500-1:1000,FC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:190 kDa; Observed MW:190 kDa

研究背景 (Background)


产品名称: TIAM2 Rabbit Monoclonal Antibody
产品货号: AMRe87589



This gene encodes a guanine nucleotide exchange factor. A highly similar mouse protein specifically activates ras-related C3 botulinum substrate 1, converting this Rho-like guanosine triphosphatase (GTPase) from a guanosine diphosphate-bound inactive state to a guanosine triphosphate-bound active state. The encoded protein may play a role in neural cell development. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of TIAM2 in K562, 293 cell lysates using TIAM2 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .