

产品名称: CIN85 Rabbit Monoclonal Antibody
产品货号: AMRe87607



产品概述 (Summary)

产品名称 (Production Name)	CIN85 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CIN85
别名 (Alternative Names)	HSB1; CIN85; GIG10; HSB-1; MIG18; CD2BP3
基因 ID (Gene ID)	30011
蛋白 ID (SwissProt ID)	Q96B97.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:73 kDa; Observed MW:82 kDa

研究背景 (Background)


产品名称: CIN85 Rabbit Monoclonal Antibody
产品货号: AMRe87607



This gene encodes an adapter protein that contains one or more N-terminal Src homology domains, a proline rich region and a C-terminal coiled-coil domain. The encoded protein facilitates protein-protein interactions and has been implicated in numerous cellular processes including apoptosis, cytoskeletal rearrangement, cell adhesion and in the regulation of clathrin-dependent endocytosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms.
[provided by RefSeq, Jul 2017]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Jurkat cells using AMRe87607 at 1:1000.

注意事项 (Note)

For research use only .