

产品名称: SSPN Rabbit Monoclonal Antibody
产品货号: AMRe87616



产品概述 (Summary)

产品名称 (Production Name)	SSPN Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SSPN
别名 (Alternative Names)	KRAG; NSPN; SPN1; SPN2; DAGA5
基因 ID (Gene ID)	8082
蛋白 ID (SwissProt ID)	Q14714.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:27 kDa; Observed MW:27 kDa

研究背景 (Background)


产品名称: SSPN Rabbit Monoclonal Antibody
产品货号: AMRe87616



This gene encodes a member of the dystrophin-glycoprotein complex (DGC). The DGC spans the sarcolemma and is comprised of dystrophin, syntrophin, alpha- and beta-dystroglycans and sarcoglycans. The DGC provides a structural link between the subsarcolemmal cytoskeleton and the extracellular matrix of muscle cells. Two alternatively spliced transcript variants that encode different protein isoforms have been described. [provided by RefSeq, Oct 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Rat heart tissue using AMRe87616 at 1:1000.

注意事项 (Note)

For research use only .