

产品名称: NMT55/p54nrb Rabbit Monoclonal Antibody
产品货号: AMRe87627



产品概述 (Summary)

产品名称 (Production Name)	NMT55/p54nrb Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NMT55/p54nrb
别名 (Alternative Names)	P54; NMT55; NRB54; MRXS34; P54NRB; PPP1R114
基因 ID (Gene ID)	4841
蛋白 ID (SwissProt ID)	Q15233.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:54 kDa; Observed MW:54 kDa

研究背景 (Background)


产品名称: NMT55/p54nrb Rabbit Monoclonal Antibody
产品货号: AMRe87627



This gene encodes an RNA-binding protein which plays various roles in the nucleus, including transcriptional regulation and RNA splicing. A rearrangement between this gene and the transcription factor E3 gene has been observed in papillary renal cell carcinoma. Alternatively spliced transcript variants have been described. Pseudogenes exist on Chromosomes 2 and 16. [provided by RefSeq, Feb 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from RAW264.7, Jurkat cells and Rat spleen tissue using AMRe87627 at 1:2000.

注意事项 (Note)

For research use only .