

产品名称: Cleaved PARP1 Rabbit Monoclonal Antibody  
产品货号: AMRe87630



## 产品概述 (Summary)

|                        |                                          |
|------------------------|------------------------------------------|
| 产品名称 (Production Name) | Cleaved PARP1 Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody   |
| 宿主 (Host)              | Rabbit                                   |
| 应用 (Application)       | WB,IHC                                   |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                          |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                             |
|------------------------|-------------------------------------------------------------|
| 基因名 (Gene Name)        | Cleaved PARP1                                               |
| 别名 (Alternative Names) | PARP; PPOL; ADPRT; ARTD1; ADPRT1; PARP-1; ADPRT 1; pADPRT-1 |
| 基因 ID (Gene ID)        | 142                                                         |
| 蛋白 ID (SwissProt ID)   | P09874.                                                     |

## 产品应用 (Application)

|                          |                                           |
|--------------------------|-------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000,IHC 1:50-1:200           |
| 蛋白分子量 (Molecular Weight) | Calculated MW:113 kDa; Observed MW:24 kDa |

## 研究背景 (Background)


**产品名称: Cleaved PARP1 Rabbit Monoclonal Antibody**  
**产品货号: AMRe87630**



This gene encodes a chromatin-associated enzyme, poly(ADP-ribosyl)transferase, which modifies various nuclear proteins by poly(ADP-ribosyl)ation. The modification is dependent on DNA and is involved in the regulation of various important cellular processes such as differentiation, proliferation, and tumor transformation and also in the regulation of the molecular events involved in the recovery of cell from DNA damage. In addition, this enzyme may be the site of mutation in Fanconi anemia, and may participate in the pathophysiology of type I diabetes. [provided by RefSeq, Jul 2008]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from Jurkat cells using AMRe87630 at 1:1000.

## 注意事项 (Note)

For research use only .