

产品名称: AF10 Rabbit Monoclonal Antibody
产品货号: AMRe87637



产品概述 (Summary)

产品名称 (Production Name)	AF10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AF10
别名 (Alternative Names)	AF10
基因 ID (Gene ID)	8028
蛋白 ID (SwissProt ID)	P55197.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,FC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:113 kDa; Observed MW:140 kDa

研究背景 (Background)


产品名称: AF10 Rabbit Monoclonal Antibody
产品货号: AMRe87637



This gene encodes a transcription factor and has been identified as a partner gene involved in several chromosomal rearrangements resulting in various leukemias. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of AF10 in HeLa, K562 cell lysates using AF10 antibody (1:1000 diluted).

注意事项 (Note)

For research use only.