

产品名称: RBPJK Rabbit Monoclonal Antibody
产品货号: AMRe87638



产品概述 (Summary)

产品名称 (Production Name)	RBPJK Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBPJK
别名 (Alternative Names)	SUH; csl; AOS3; CBF1; KBF2; RBP-J; RBPJK; IGKJRB; RBPSUH; IGKJRB1
基因 ID (Gene ID)	3516
蛋白 ID (SwissProt ID)	Q06330.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:100-1:200, FC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:56 kDa; Observed MW:61 kDa

研究背景 (Background)


产品名称: RBPJK Rabbit Monoclonal Antibody
产品货号: AMRe87638




The protein encoded by this gene is a transcriptional regulator important in the Notch signaling pathway. The encoded protein acts as a repressor when not bound to Notch proteins and an activator when bound to Notch proteins. It is thought to function by recruiting chromatin remodeling complexes containing histone deacetylase or histone acetylase proteins to Notch signaling pathway genes. Several transcript variants encoding different isoforms have been found for this gene, and several pseudogenes of this gene exist on chromosome 9. [provided by RefSeq, Oct 2013]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of RBPJK in Rat Brain, C6 cell lysates using RBPJK antibody (1:1000 diluted).



Western blot analysis of extracts from RAW264.7 cells using AMRe87638 at 1:3000.

注意事项 (Note)

For research use only.