

产品名称: CD3 Rabbit Monoclonal Antibody
产品货号: AMRe87639



产品概述 (Summary)

产品名称 (Production Name)	CD3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD3
别名 (Alternative Names)	T3
基因 ID (Gene ID)	915, 916, 917, 919
蛋白 ID (SwissProt ID)	P04234, P07766, P09693, P20963 .

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:200,FC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:23 kDa; Observed MW:23 kDa

研究背景 (Background)


产品名称: CD3 Rabbit Monoclonal Antibody
产品货号: AMRe87639



CD3 is a complex of proteins that directly associates with the T cell receptor (TCR). The TCR/CD3 complex of T-lymphocytes consists of either a TCR alpha/beta or TCR gamma/delta heterodimer coexpressed at the cell surface with the invariant subunits of CD3 labeled gamma, delta, epsilon, zeta, and eta. The TCR recognizes antigens bound to major histocompatibility complex (MHC) molecules. TCR-mediated peptide-MHC recognition is transmitted to the CD3 complex, leading to the intracellular signal transduction. CD3 is considered to be a pan-T cell marker for detection of normal and neoplastic T cells.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Jurkat cells using AMRe87639 at 1:1000.

注意事项 (Note)

For research use only .