

产品名称: USP16 Rabbit Monoclonal Antibody
产品货号: AMRe87641



产品概述 (Summary)

产品名称 (Production Name)	USP16 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP16
别名 (Alternative Names)	UBPM; 1200004E02Rik; 2810483I07Rik; 6330514E22Rik
基因 ID (Gene ID)	74112, 288306
蛋白 ID (SwissProt ID)	Q99LG0, Q2KJ09.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:93 kDa; Observed MW:110 kDa

研究背景 (Background)


产品名称: USP16 Rabbit Monoclonal Antibody
产品货号: AMRe87641



Predicted to enable several functions, including cysteine-type peptidase activity; ubiquitin binding activity; and zinc ion binding activity. Involved in monoubiquitinated histone H2A deubiquitination and positive regulation of transcription by RNA polymerase II. Predicted to be located in nucleus. Is expressed in several structures, including alimentary system; genitourinary system; integumental system; nervous system; and sensory organ. Orthologous to human USP16 (ubiquitin specific peptidase 16). [provided by Alliance of Genome Resources, Apr 2022]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse testis tissue using AMRe87641 at 1:1000.

注意事项 (Note)

For research use only .