

产品名称: TYRP1 Rabbit Monoclonal Antibody
产品货号: AMRe87645



产品概述 (Summary)

产品名称 (Production Name)	TYRP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TYRP1
别名 (Alternative Names)	TRP; CAS2; CATB; GP75; OCA3; TRP1; TYRP; b-PROTEIN
基因 ID (Gene ID)	7306, 22178
蛋白 ID (SwissProt ID)	P17643, P07147.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:200-1:1000,ICC/IF 1:100-1:200,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:61 kDa; Observed MW:75 kDa

研究背景 (Background)


产品名称: TYRP1 Rabbit Monoclonal Antibody
产品货号: AMRe87645



This gene encodes a melanosomal enzyme that belongs to the tyrosinase family and plays an important role in the melanin biosynthetic pathway. Defects in this gene are the cause of rufous oculocutaneous albinism and oculocutaneous albinism type III. [provided by RefSeq, Mar 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from B16-F10 cells using AMRe87645 at 1:5000.

注意事项 (Note)

For research use only .