

产品名称: Midkine Rabbit Monoclonal Antibody
产品货号: AMRe87661



产品概述 (Summary)

产品名称 (Production Name)	Midkine Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Midkine
别名 (Alternative Names)	MK; ARAP; NEGF2
基因 ID (Gene ID)	4192
蛋白 ID (SwissProt ID)	P21741.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,ICC/IF 1:100-1:200,FC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:16 kDa; Observed MW:16 kDa

产品名称: Midkine Rabbit Monoclonal Antibody
产品货号: AMRe87661




研究背景 (Background)


This gene encodes a member of a small family of secreted growth factors that binds heparin and responds to retinoic acid. The encoded protein promotes cell growth, migration, and angiogenesis, in particular during tumorigenesis. This gene has been targeted as a therapeutic for a variety of different disorders. Alternatively spliced transcript variants encoding multiple isoforms have been observed. [provided by RefSeq, Jul 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Midkine in U251,293 cell lysates using Midkine antibody(1:1000 diluted).



Western blot analysis of extracts from SH-SY5Y cells using AMRe87661 at 1:1000.

注意事项 (Note)

For research use only .