

产品名称: Inhibin beta A Rabbit Monoclonal Antibody
产品货号: AMRe87664



产品概述 (Summary)

产品名称 (Production Name)	Inhibin beta A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Inhibin beta A
别名 (Alternative Names)	EDF; FRP
基因 ID (Gene ID)	3624
蛋白 ID (SwissProt ID)	P08476.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:47 kDa; Observed MW:44 kDa

研究背景 (Background)


产品名称: Inhibin beta A Rabbit Monoclonal Antibody
产品货号: AMRe87664




This gene encodes a member of the TGF-beta (transforming growth factor-beta) superfamily of proteins. The encoded preproprotein is proteolytically processed to generate a subunit of the dimeric activin and inhibin protein complexes. These complexes activate and inhibit, respectively, follicle stimulating hormone secretion from the pituitary gland. The encoded protein also plays a role in eye, tooth and testis development. Elevated expression of this gene may be associated with cancer cachexia in human patients. [provided by RefSeq, Aug 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Inhibin beta A in Rat Brain, C6, CHO-K1, HeLa cell lysates using Inhibin beta A antibody (1:1000 diluted).



Western blot analysis of extracts from RAW264.7 cells using AMRe87664 at 1:3000.

注意事项 (Note)

For research use only .