

产品名称: PPP2CA/B Rabbit Monoclonal Antibody
产品货号: AMRe87673



产品概述 (Summary)

产品名称 (Production Name)	PPP2CA/B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPP2CA/B
别名 (Alternative Names)	RP-C; PP2Ac; PP2CA; PP2Calpha
基因 ID (Gene ID)	5515, 5516
蛋白 ID (SwissProt ID)	P62714, P67775.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:36 kDa; Observed MW:36 kDa

研究背景 (Background)


产品名称: PPP2CA/B Rabbit Monoclonal Antibody
产品货号: AMRe87673




This gene encodes the phosphatase 2A catalytic subunit. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. This gene encodes an alpha isoform of the catalytic subunit. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Jurkat cells using AMRe87673 at 1:1000.



Western blot analysis of extracts from HepG2, C2C12 cells and Rat spleen tissue using AMRe87673 at 1:1000.

注意事项 (Note)

For research use only .