

产品名称: ANP32B Rabbit Monoclonal Antibody
产品货号: AMRe87678



产品概述 (Summary)

产品名称 (Production Name)	ANP32B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANP32B
别名 (Alternative Names)	APRIL; SSP29; PHAPI2
基因 ID (Gene ID)	10541, 67628, 170724
蛋白 ID (SwissProt ID)	Q92688, Q9EST5, Q9EST6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:200-1:500,FC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:29 kDa; Observed MW:29 kDa

研究背景 (Background)


产品名称: ANP32B Rabbit Monoclonal Antibody
产品货号: AMRe87678



Enables RNA polymerase binding activity and histone binding activity. Involved in several processes, including activation of cysteine-type endopeptidase activity involved in apoptotic process; nucleosome assembly; and positive regulation of protein export from nucleus. Located in cytoplasm and nucleoplasm. Colocalizes with nucleolus. [provided by Alliance of Genome Resources, Apr 2022]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Jurkat cells using AMRe87678 at 1:1000.

注意事项 (Note)

For research use only .