

产品名称: PCSK1 Rabbit Monoclonal Antibody
产品货号: AMRe87681



产品概述 (Summary)

产品名称 (Production Name)	PCSK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PCSK1
别名 (Alternative Names)	PC1; PC3; NEC1; SPC3; PC1/3; BMIQ12
基因 ID (Gene ID)	5122, 18548, 25204
蛋白 ID (SwissProt ID)	P29120, P63239, P28840.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:100-1:200,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:84 kDa; Observed MW:70 kDa

研究背景 (Background)


产品名称: PCSK1 Rabbit Monoclonal Antibody
产品货号: AMRe87681




This gene encodes a member of the subtilisin-like proprotein convertase family, which includes proteases that process protein and peptide precursors trafficking through regulated or constitutive branches of the secretory pathway. The encoded protein undergoes an initial autocatalytic processing event in the ER to generate a heterodimer which exits the ER and sorts to subcellular compartments where a second autocatalytic even takes place and the catalytic activity is acquired. The protease is packaged into and activated in dense core secretory granules and expressed in the neuroendocrine system and brain. This gene encodes one of the seven basic amino acid-specific members which cleave their substrates at single or paired basic residues. It functions in the proteolytic activation of polypeptide hormones and neuropeptides precursors. Mutations in this gene have been associated with susceptibility to obesity and proprotein convertase 1/3 deficiency. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene [provided by RefSeq, Jan 2014]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse small intestine tissue using AMRe87681 at 1:1000.



Western blot analysis of extracts from Mouse brain tissue using AMRe87681 at 1:100.

产品名称: PCSK1 Rabbit Monoclonal Antibody
产品货号: AMRe87681



注意事项 (Note)

For research use only .