

产品名称: Phospho-AP2M1 (Thr156) Rabbit Monoclonal Antibody
产品货号: AMRe87683



产品概述 (Summary)

产品名称 (Production Name)	Phospho-AP2M1 (Thr156) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-AP2M1 (Thr156)
别名 (Alternative Names)	mu2; AP50; CLAPM1
基因 ID (Gene ID)	1173
蛋白 ID (SwissProt ID)	Q96CW1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 50 kDa

产品名称: Phospho-AP2M1 (Thr156) Rabbit Monoclonal Antibody
产品货号: AMRe87683




研究背景 (Background)


This gene encodes a subunit of the heterotetrameric coat assembly protein complex 2 (AP2), which belongs to the adaptor complexes medium subunits family. The encoded protein is required for the activity of a vacuolar ATPase, which is responsible for proton pumping occurring in the acidification of endosomes and lysosomes. The encoded protein may also play an important role in regulating the intracellular trafficking and function of CTLA-4 protein. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Phospho-AP2M1 (Thr156) in K562 cell lysates using Phospho-AP2M1 (Thr156) antibody(1:1000 diluted).



Western blot analysis of extracts from HeLa, PC-12 cells using AMRe87683 at 1:1000.

注意事项 (Note)

For research use only .