

产品名称: Transcription Initiation Factor IIB Rabbit
Monoclonal Antibody
产品货号: AMRe87684



产品概述 (Summary)

产品名称 (Production Name)	Transcription Initiation Factor IIB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Transcription Initiation Factor IIB
别名 (Alternative Names)	TF2B; TFIIB
基因 ID (Gene ID)	2959
蛋白 ID (SwissProt ID)	Q00403.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:35 kDa; Observed MW:35 kDa


**产品名称: Transcription Initiation Factor IIB Rabbit
Monoclonal Antibody**
产品货号: AMRe87684

研究背景 (Background)

This gene encodes the general transcription factor IIB, one of the ubiquitous factors required for transcription initiation by RNA polymerase II. The protein localizes to the nucleus where it forms a complex (the DAB complex) with transcription factors IID and IIA. Transcription factor IIB serves as a bridge between IID, the factor which initially recognizes the promoter sequence, and RNA polymerase II. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse muscle tissue using AMRe87684 at 1:1000.

注意事项 (Note)

For research use only .