

产品名称: BTG1 Rabbit Monoclonal Antibody
产品货号: AMRe87710



产品概述 (Summary)

产品名称 (Production Name)	BTG1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BTG1
别名 (Alternative Names)	APRO2
基因 ID (Gene ID)	694, 12226, 29618
蛋白 ID (SwissProt ID)	P62324, P62325, Q63073.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:19 kDa; Observed MW:19 kDa

研究背景 (Background)


产品名称: BTG1 Rabbit Monoclonal Antibody
产品货号: AMRe87710



This gene is a member of an anti-proliferative gene family that regulates cell growth and differentiation. Expression of this gene is highest in the G0/G1 phases of the cell cycle and downregulated when cells progressed through G1. The encoded protein interacts with several nuclear receptors, and functions as a coactivator of cell differentiation. This locus has been shown to be involved in a t(8;12)(q24;q22) chromosomal translocation in a case of B-cell chronic lymphocytic leukemia. [provided by RefSeq, Oct 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from SH-SY5Y cells using AMRe87710 at 1:1000.

注意事项 (Note)

For research use only .