

产品名称: BRN3A Rabbit Monoclonal Antibody
产品货号: AMRe87711



产品概述 (Summary)

产品名称 (Production Name)	BRN3A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BRN3A
别名 (Alternative Names)	BRN3A; RDC-1; Oct-T1; brn-3A
基因 ID (Gene ID)	5457, 18996, 114503
蛋白 ID (SwissProt ID)	Q01851, P17208, P20266.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:43 kDa; Observed MW:43 kDa

研究背景 (Background)


产品名称: BRN3A Rabbit Monoclonal Antibody
产品货号: AMRe87711



This gene encodes a member of the POU-IV class of neural transcription factors. This protein is expressed in a subset of retinal ganglion cells and may be involved in the developing sensory nervous system. This protein may also promote the growth of cervical tumors. A translocation of this gene is associated with some adult acute myeloid leukemias. [provided by RefSeq, Mar 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from T24 cells using AMRe87711 at 1:1000.

注意事项 (Note)

For research use only .