

产品名称: Bovine Serum Albumin Rabbit Monoclonal Antibody
产品货号: AMRe87712



产品概述 (Summary)

产品名称 (Production Name)	Bovine Serum Albumin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Bovine Serum Albumin
别名 (Alternative Names)	HSA; FDAHT; PRO0883; PRO0903; PRO1341
基因 ID (Gene ID)	280717, 213, 11657, 24186
蛋白 ID (SwissProt ID)	P02769, P02768, P07724, P02770.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:69 kDa; Observed MW:69 kDa

产品名称: Bovine Serum Albumin Rabbit Monoclonal Antibody
产品货号: AMRe87712




研究背景 (Background)

This gene encodes the most abundant protein in human blood. This protein functions in the regulation of blood plasma colloid osmotic pressure and acts as a carrier protein for a wide range of endogenous molecules including hormones, fatty acids, and metabolites, as well as exogenous drugs. Additionally, this protein exhibits an esterase-like activity with broad substrate specificity. The encoded preproprotein is proteolytically processed to generate the mature protein. A peptide derived from this protein, EPI-X4, is an endogenous inhibitor of the CXCR4 chemokine receptor. [provided by RefSeq, Jul 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe87712 at 1:1000.

注意事项 (Note)

For research use only .