

产品名称: NIT1 Rabbit Monoclonal Antibody
产品货号: AMRe87714



产品概述 (Summary)

产品名称 (Production Name)	NIT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NIT1
别名 (Alternative Names)	Deaminated glutathione amidase; dGSH amidase; Nitrilase homolog 1
基因 ID (Gene ID)	4817
蛋白 ID (SwissProt ID)	Q86X76.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:20-1:50,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:36 kDa; Observed MW:35 kDa

研究背景 (Background)


产品名称: NIT1 Rabbit Monoclonal Antibody
产品货号: AMRe87714



This gene encodes a member of the nitrilase protein family with homology to bacterial and plant nitrilases, enzymes that cleave nitriles and organic amides to the corresponding carboxylic acids plus ammonia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of NIT1 in HeLa,293T,Jurkat,3T3 cell lysates using NIT1 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .