

产品名称: ARHGAP25 Rabbit Monoclonal Antibody
产品货号: AMRe87716



产品概述 (Summary)

产品名称 (Production Name)	ARHGAP25 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARHGAP25
别名 (Alternative Names)	KAIA0053; HEL-S-308
基因 ID (Gene ID)	9938
蛋白 ID (SwissProt ID)	P42331.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:73 kDa; Observed MW:73 kDa

研究背景 (Background)


产品名称: ARHGAP25 Rabbit Monoclonal Antibody
产品货号: AMRe87716



ARHGAPs, such as ARHGAP25, encode negative regulators of Rho GTPases (see ARHA; MIM 165390), which are implicated in actin remodeling, cell polarity, and cell migration (Katoh and Katoh, 2004 [PubMed 15254788]).[supplied by OMIM, Mar 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of ARHGAP25 in Rat brain, Ramos, Jurkat cell lysates using ARHGAP25 antibody (1:1000 diluted).

注意事项 (Note)

For research use only.