

产品名称: SMARCC2/BAF170 Rabbit Monoclonal Antibody
产品货号: AMRe87719



产品概述 (Summary)

产品名称 (Production Name)	SMARCC2/BAF170 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMARCC2/BAF170
别名 (Alternative Names)	CSS8; Rsc8; BAF170; CRACC2
基因 ID (Gene ID)	6601, 68094, 362815
蛋白 ID (SwissProt ID)	Q8TAQ2, Q6PDG5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC/IF 1:100-1:200,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:133 kDa; Observed MW:170 kDa

产品名称: SMARCC2/BAF170 Rabbit Monoclonal Antibody
产品货号: AMRe87719




研究背景 (Background)


The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI and contains a predicted leucine zipper motif typical of many transcription factors. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from A549 cells using AMRe87719 at 1:1000.



Western blot analysis of extracts from HCT116 cells and Rat spleen tissue using AMRe87719 at 1:1000.

注意事项 (Note)

For research use only .