

产品名称: C3 Rabbit Monoclonal Antibody
产品货号: AMRe87723



产品概述 (Summary)

产品名称 (Production Name)	C3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C3
别名 (Alternative Names)	ASP; C3a; C3b; AHUS5; ARMD9; CPAMD1; HEL-S-62p
基因 ID (Gene ID)	718
蛋白 ID (SwissProt ID)	P01024.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:187 kDa; Observed MW:187 kDa

研究背景 (Background)


产品名称: C3 Rabbit Monoclonal Antibody
产品货号: AMRe87723



Complement component C3 plays a central role in the activation of complement system. Its activation is required for both classical and alternative complement activation pathways. The encoded preproprotein is proteolytically processed to generate alpha and beta subunits that form the mature protein, which is then further processed to generate numerous peptide products. The C3a peptide, also known as the C3a anaphylatoxin, modulates inflammation and possesses antimicrobial activity. Mutations in this gene are associated with atypical hemolytic uremic syndrome and age-related macular degeneration in human patients. [provided by RefSeq, Nov 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from T24 cells using AMRe87723 at 1:1000.

注意事项 (Note)

For research use only .