

产品名称: ALDH16A1 Rabbit Monoclonal Antibody
产品货号: AMRe87728



产品概述 (Summary)

产品名称 (Production Name)	ALDH16A1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALDH16A1
别名 (Alternative Names)	ALDH16A1
基因 ID (Gene ID)	126133, 69748, 361571
蛋白 ID (SwissProt ID)	Q8IZ83, Q571I9, Q3T1L0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:85 kDa; Observed MW:85 kDa

研究背景 (Background)


产品名称: ALDH16A1 Rabbit Monoclonal Antibody
产品货号: AMRe87728



This gene encodes a member of the aldehyde dehydrogenase superfamily. The family members act on aldehyde substrates and use nicotinamide adenine dinucleotide phosphate (NADP) as a cofactor. This gene is conserved in chimpanzee, dog, cow, mouse, rat, and zebrafish. The protein encoded by this gene interacts with maspardin, a protein that when truncated is responsible for Mast syndrome. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2, HeLa cells using AMRe87728 at 1:1000.

注意事项 (Note)

For research use only .