

产品名称: AATF Rabbit Monoclonal Antibody
产品货号: AMRe87732



产品概述 (Summary)

产品名称 (Production Name)	AATF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AATF
别名 (Alternative Names)	DED; BFR2; CHE1; CHE-1
基因 ID (Gene ID)	26574, 56321
蛋白 ID (SwissProt ID)	Q9NY61, Q9JKX4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:63 kDa; Observed MW:80 kDa

研究背景 (Background)


产品名称: AATF Rabbit Monoclonal Antibody
产品货号: AMRe87732




The protein encoded by this gene was identified on the basis of its interaction with MAP3K12/DLK, a protein kinase known to be involved in the induction of cell apoptosis. This gene product contains a leucine zipper, which is a characteristic motif of transcription factors, and was shown to exhibit strong transactivation activity when fused to Gal4 DNA binding domain. Overexpression of this gene interfered with MAP3K12 induced apoptosis. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87732 at 1:1000.



Western blot analysis of extracts from HepG2, A431 cells using AMRe87732 at 1:1000.

注意事项 (Note)

For research use only .